


The Kronner Manipujector®

Intra-operatively Adjustable
Uterine Manipulator/Injector

Features and Benefits

- Affords rapid intra-operative repositioning
- Provides excellent manipulative control
- Securely seals the internal os for a no-leak injection
- Safe, simple, and economical



The Kronner Manipujector[®] is the standard in disposable, balloon-tipped uterine manipulator/injectors. Designed for patient safety, ease-of-use, convenience, and maximum procedural versatility.

Unique Pinch-grip Handle for Rapid, Effortless Intra-Operative Adjustment

The Kronner Manipujector's unique glide action, pinch-grip handle is designed to eliminate the time-consuming tasks required to remove, reinsert, and reinflate when intra-operative repositioning is required. Just pinch the grip and glide the handle to the new position. When in place, simply release the grip, which sets the handle firmly, and continue with the procedure.

Anatomically Curved Shape Provides Excellent Manipulative Control

The radial design optimizes cervical visualization, allowing easy insertion and maneuverability of either anteverted or retroverted uteri. The ergonomic curvature affords superb elevation for ease of surgical dissection and visualization.

Positive Balloon Positioning Reduces Potential for Expulsion

The elimination of spring-controlled automatic retraction obviates the need to maintain forward pressure during balloon inflation and helps avoid inadvertent balloon inflation within the cervical canal – a common cause of spontaneous balloon expulsion requiring device reinsertion.

Safe, simple, and economical to use, the Kronner Manipujector is designed to afford maximum operator convenience and minimal uterine trauma.

Ordering Information

Code	Description
Kronner Manipujector	
6003	Kronner Manipujector [®] (12 per box)



Proudly distributed by



Endotherapeutics
IMPROVING HEALTHCARE



1300 00 ENDO (3636)
info@endotherapeutics.com.au
www.endotherapeutics.com.au
Suite 301, 2 Banfield Road
Macquarie Park NSW 2113



CooperSurgical is the leading company dedicated to providing medical devices and procedure solutions that improve health care delivery to women regardless of clinical setting. Our company is fostering that position through expansion of its core businesses and introduction of advanced technology-based products, which aid clinicians in the management and treatment of commonly seen conditions.

800.243.2974 | 203.601.5200 | coopersurgical.com | 95 Corporate Drive | Trumbull, CT 06611
Also on [YouTube.com/coopersurgical](https://www.youtube.com/coopersurgical)