



# Lone Star APS Retractor™

For Obstetric and Gynecologic Procedures

See the difference unobstructed  
transvaginal surgery can make

- Unobstructed visualization and access
- Designed specifically for modern pelvic and vaginal surgery
- Effectively manages the surgical site, simplifying the procedure
- VELCRO® backing for increased stability and ease of placement
- Frees up OR personnel to perform other tasks





Convenient clip manages the Foley catheter

VELCRO backing facilitates site setup and management

Closely placed slots allow greater precision in hook/stay placement

Custom holes for non-perforating towel clamps keep retractor stable and stationary

Ideal width for a variety of instrument introduction

Unique, open-end design facilitates use of weighted speculum, or access to the posterior compartment

### The APS Retractor designed specifically for:

- Reconstructive pelvic surgery
  - Sling procedures
  - Anterior repair procedures
  - Posterior repair procedures
- Mesh excision
- Episiotomy repair
- Vaginal hysterectomy
- Vaginal biopsies

## Lone Star APS Retractor™ For Obstetric and Gynecologic Procedures

The APS Retractor combined with the recommended Stays provides excellent visual exposure of the vagina.

### Unprecedented access – for surgical excellence

The versatile APS system provides unparalleled exposure. For simpler, better-managed surgeries, the innovative design of the self-retaining APS Retractor™ offers a clear pathway, unobstructed anatomical visibility and dramatic improvement in site management with minimum effort. Adjustable clips, stable VELCRO® backing and ultra-precise hook placements quickly adapt the APS Retractor to present the surgical site exactly as needed. The system frees up OR personnel, saving time and operating expense. There's nothing like it – and nothing that performs as well.

## Lone Star APS Retractor™ Ordering Information

Code	Description
<b>APS Retractor (Packaged singly with 2 catheter clips)</b>	
3715	Height 16.6 cm x Width 16.6 cm, Weight 77g
<b>Recommended Stays</b>	
3311-8G	5 mm Sharp Hook, 8/Pack
3350-8G	12 mm Blunt Hook, 8/Pack
3327-4G	7 mm Two-finger Sharp Hook, 4/Pack
3384-4	12 mm Two-finger Blunt Hook, 4/Pack

Products are sterile, single-use only



Endotherapeutics

### Distributed by

**Endotherapeutics**

ABN 89 088 045 103

Level 1, 51 Rawson St

Epping NSW 2121 Australia

T: 1300 00 3636 F: +61 (2) 9869 7723

E: [info@endotherapeutics.com.au](mailto:info@endotherapeutics.com.au)

[www.endotherapeutics.com.au](http://www.endotherapeutics.com.au)



CooperSurgical is the leading company dedicated to providing medical devices and procedure solutions that improve health care delivery to women regardless of clinical setting. Our company is fostering that position through expansion of its core businesses and introduction of advanced technology-based products, which aid clinicians in the management and treatment of commonly seen conditions.

To find out more about the Lone Star APS Retractor or to place an order Call 800.243.2974 or visit [www.coopersurgical.com](http://www.coopersurgical.com).

800.243.2974 | 203.601.5200 | [coopersurgical.com](http://coopersurgical.com) | 95 Corporate Drive | Trumbull, CT 06611  
Also on [YouTube.com/coopersurgical](https://www.youtube.com/coopersurgical)